USER GUIDE 1-7

GUIA DEL USUARIO 8-14

COMODITÀ

Excellence in design and safety

COMODITÀ Siena
all terrain

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model number: COM 700

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Thank you for choosing the COMODITÀ Siena All Terrain Rolling Walker (Ref. COM 700). At COMODITÀ, we are committed to your comfort, safety, and mobility. Our products are designed to make daily activities easier and to help you move with confidence and independence.

The Siena All Terrain Rolling Walker combines a modern, patented design with strong, reliable construction. Tested and approved under international safety standards, it supports users up to 250 lbs. and is built with high-quality materials to provide the stability and durability you can trust.

This model includes a removable cup holder, cane holder, and washable nylon shopping bag, adding convenience to your everyday use. With two 12" front wheels and two 10" rear wheels, made of soft, solid PU tires, the Siena All Terrain is ideal for both indoor and outdoor use in all kinds of terrains.

Please do not hesitate to contact us at info@comoditausea.com if you have any doubts, questions or problems. We are here to help!

Main Features:

- Extra-large 12" front wheels and 10" rear wheels with solid, soft PU tires, for safer indoor and outdoor use
- Wide, comfortable nylon seat
- Adjustable handle height
- Padded ergonomic reinforced backrest, with built-in phone pocket
- Convenient removable cup holder
- Cane holder and included
- Hand brake locks for safety while sitting
- Sideways folding construction makes it very compact and easy to transport
- Two light reflectors at the front and two at the back for additional safety
- Nylon shopping bag with light reflecting strap for added security

Product Dimensions and Specifications

- Unit weight: 24 lb. (11 kg)
- Weight capacity: 250 lb. (114 Kg)
- Recommended user height: 5'4" to 6'4" (1,62 m to 1,90 m)
- Overall height: 36" / 44" (91 cm / 111 cm)
- Overall length: 31" (79 cm)
- Overall width: 26" (66 cm)
- Folded measurements: 36"x 12"x 31"
- Front wheel diameter: 12" (30 cm)
- Back wheel diameter: 10" (20 cm)
- Seat height from floor: 23" (59 cm)
- Seat depth: 11" (28 cm)
- Seat width: 18" (46 cm)
- Size of nylon bag: 15" W, 8" D, 10" H (38 cm W, 20cm D, 25cm H)

Assembly Instructions

Assembling your COMODITÀ Siena All Terrain Rolling Walker is extremely simple!

1. Open the walker and lock it in the opened position by pushing down the sides of the seat onto the main frame until you hear a click.

2. Insert both sides of the nylon shopping bag straps around the nylon bag supports on the sides of the seat.

3. Insert the cup holder into its plastic support on the side of the walker.

4. Adjust the height of the handles by standing upright behind the walker with the shoulders relaxed and arms hanging at your side. The handgrips should be at a height even with your wrists.

5. Press the orange button to release the handle and set it to your desired height.

IT IS THAT SIMPLE! NO TOOLS, NO LOOSE PARTS!

Operating Instructions

TO WALK:

- Position yourself behind the walker and hold both handles firmly.
- Roll your walker slowly forward so that the rear wheels are a few inches ahead of your body.
- While maintaining a steady and safe roll forward, place your foot in line with the rear wheels.
- Repeat the movement with your other foot.
- To slow down, slowly squeeze the brake handles up.
- To lock the brakes, push down on the brake handles until the brakes "click" into locked position. Do not push it any further than this point.
- To unlock the brakes and resume walking, pull up the brake handles until you hear a "click."

TO USE THE SEAT OF YOUR WALKER:

- Make sure your walker is in a flat surface, not on an incline.
- LOCK the brakes before sitting.
- Make sure your feet rest flat on the ground.

Care and Maintenance

PERIODIC BRAKE ADJUSTMENT

To adjust the sensitivity of the brake:

- To tighten the brake (reduce the gap between the brake pad and wheel) rotate the adjusting nut clockwise.
- To loosen the brake (increase the gap between the brake and the wheel) rotate the adjusting nut counterclockwise.

WARNING!

- Over adjustment of the brake might prevent the wheels from moving freely.
- Under adjustment of the brake may cause the brake pad not to engage at all.
- **DO NOT USE** THE WALKER IF THE BRAKES ARE NOT WORKING PROPERLY.

Safety Instructions

WARNING!!! Failure to follow safety instructions may result in serious injuries.

- **DO NOT** use this product as a wheelchair.
- **DO NOT** transport someone while sitting on this walker.
- **DO NOT** exceed maximum weight capacity.
- Brakes **MUST** be on **LOCKED** position before sitting on the seat.
- **DO NOT** self-propel the walker while seated.
- **DO NOT** sit on the walker on an incline.
- **INSPECT BRAKES** before each use.
- **DO NOT** use the walker if the brakes are not working properly.
- **ALL WHEELS MUST** be in contact with the floor or ground at all times.
- **BACKREST IS NOT** intended to support your full weight. It is intended to provide back support while you are seated.
- **DO NOT** use in escalators, stairs, moving sidewalks or uneven surfaces.
- **BOTH FEET** should rest flat on the ground when sitting.
- **ALWAYS** check if the seat is locked down and well connected to the frame when the rollator is open.

Warranty

Your Comodità Rolling Walker is warranted to be free of defects in material and workmanship for a period of one year from the date of purchase for the original owner. This warranty does not cover product failure due to user negligence, misuse of the product, failure to adhere to instructions, or normal wear and tear. In the event of a defect covered by this warranty, we will, at our discretion, repair or replace the product. This warranty does not include any labor or shipping charges. Please contact our customer service at info@comoditausa.com. We will provide a return authorization number. You will be required to include a copy of your original bill of sale.

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COMODITÀ LLC USA