

DISTINCT PRODUCTS



EXPERIENCES THAT MATTER



PAEDIATRIC AND NEONATAL SECUREMENT

High Flow Cannulaide® Securement Device	C2
CPAP Cannulaide® Securement Device.....	C4
Cannula and Tubing Securement	C6
General Purpose Securement.....	C7
Tilson Trach Guard.....	C8

HIGH FLOW CANNULAIDE® SECUREMENT DEVICE

Protect the skin, keep the leak. Designed for your tiniest babies who need high flow therapy, with a gentle touch. One system combines both securement and protection. Made with soft silicone material to avoid skin irritation and align with NICU skin care guidelines. Inspired by seeing the challenges that NICU clinicians face at the bedside.




HIGH FLOW CANNULAIDE® SECUREMENT DEVICE



High Flow Cannulaide® Securement Device for Paediatric and Neonatal Use

Protects the septum and nares to maintain contact between skin and cannula. Clear, see-through material allows for easier visualization of skin integrity and repositioning, if needed. Two sizes for high flow nasal cannula use to reduce worry of a seal or occlusion.



Description	Units per case	Item number
High Flow Cannulaide silicone cannula securement device for infants <1250 g (0.84" x 3.00")	Box of 25	HF-101-0-25
High Flow Cannulaide silicone cannula securement device for infants >1250 g (0.813" x 3.50")	Box of 25	HF-102-0-25

High Flow Cannulaide securement device is not for use with non-invasive ventilation (e.g. CPAP, bubble CPAP, etc.).

CPAP CANNULAIDE® SECUREMENT DEVICE

Protect the skin, enhance the seal. Gentle hydrocolloid material provides a unique protective barrier between the cannula and baby's delicate skin, providing a more consistent seal, while reducing the need to "upsized" prongs due to nasal leak.



CPAP CANNULAIDE® SECUREMENT DEVICE

CPAP Cannulaide® Securement Device For Paediatric and Neonatal Use

- Soft, low-profile hook and loop strips provide for easy securement and repositioning of the nasal interface
- Hook and loop system prevents the cannula from disconnecting from patient movement³
- Works well in humid conditions, even up to 24 hours of continuous use³



Description	Units per case	Item number
CPAP Cannulaide® Securement Device, size 0, for infants <700 g	Box of 25	CA100-0-25
CPAP Cannulaide® Securement Device, size 1, for infants 700-1250 g	Box of 25	CA101-0-25
CPAP Cannulaide® Securement Device, size 2, for infants 1251-2000 g	Box of 25	CA102-0-25
CPAP Cannulaide® Securement Device, size 3, for infants 2001-3000 g	Box of 25	CA103-0-25
CPAP Cannulaide® Securement Device, size 4, for infants >3000 g	Box of 25	CA104-0-25
CPAP Cannulaide® Securement Device, Variety Pack, 3 qty, size 0; 7 qty, size 1; 10 qty, size 2; 3 qty, size 3; 2 qty, size 4	Box of 25	CA110-0-1

CPAP Cannulaide® Securement Device is not for use with low- or high-flow oxygen cannulas, nCPAP masks, or endotracheal tubes.

CANNULA AND TUBING SECUREMENT

Mini-Whiskers™ Securement Device

Ideal for premature infants less than 2000 g. A gentle securement system for nasal cannulas and feeding tubes that uses hydrocolloid material designed to provide gentle adhesion for 24 to 36 hours. Velcro strips help keep tubing in position, reduce skin friction, and allow for easy adjustments when needed.




Sticky Whiskers™ Securement Device

A gentle securement system for nasal cannulas and feeding tubes, appropriate for neonate to paediatric patients. Made with hydrocolloid material designed to provide gentle adhesion for 24 to 36 hours. Uses Velcro strips to help keep tubing in position, reduce skin friction, and allow for easy adjustments when needed.



Tender Grips[®]

Useful for infants to adult patients. Soft, skin-toned acrylic material designed to hold cannula tubing in position, generally used around the patient's cheeks. An internal flap allows cannula adjustment without removing the securement from the patient's skin.



Description	Units per case	Item number
Mini-Whiskers Securement Device exclusively for early-premature infants (<2000 g) 6 - 0.5" x 1.4" per pouch, including 12 securement loops	20 pouches per box	MW100-0-120
Sticky Whiskers Securement Device one size fits all for infants 1 - 0.436" x 3.5" per pouch, including 2 securement loops	Box of 25	SW102-0-25
Tender Grips skin fixation system	25 pairs	1005-0-25
Tender Grips skin fixation system	100 pairs	1006-0-100

GENERAL PURPOSE SECUREMENT

Hydrocolloid Strips

For all patients. Pre-cut and ready to apply. Hydrocolloid material provides gentle adhesion to the skin for 24 to 36 hours. Also ideal for securing feeding tubes, EKG wires, and other loose tubing.



1217-0-25

Silicone Strips

For all patients, especially premature. Made with a delicate silicone material that resists moisture and can be easily repositioned and cleanly detached from the skin without causing skin irritations. Also ideal for securing feeding tubes, EKG wires, and other loose tubing.



SS100

Grippies Fastener

A non-adhesive hook and loop tape for bundling tubing and securing it to other devices. Perfect for organising cords and tubing around a patient.



GR101

Description	Units per case	Item number
Hydrocolloid strips with 2 clear securement strips per sheet	25 sheets	1217-0-25
All-purpose, moisture resistant Silicone Strips (8 per sheet, 5.25" x 6.00")	10 sheets per box	SS100P-0-80
All-purpose, moisture resistant Silicone Strips (8 per sheet, 5.25" x 6.00"), with 10 sheets per box	20 boxes per case	SS100C-0-1600
Grippies non-toxic hook and loop fastener strip (9' roll x 0.50" wide)	5 rolls	GR101P-0-5
Grippies non-toxic hook and loop fastener strip (9' roll x 0.50" wide)	150 rolls	GR101C-0-150

Our most popular models are indicated by **Salter Green** throughout this catalog.

TILSON TRACH GUARD

Tracheostomy Tube Protection

For paediatric patients, it protects an exposed tracheal tube from accidental occlusion by helping to prevent blockage of the airway from bedding, dressings, clothing, or chubby baby chins. Fits securely to standard 15 mm tube connectors, with a large central hole that allows for ample suctioning. Bright blue for easy identification.



Description

Units

Item number

Tilson Trach Guard

25 per box

TP100-0-25

For more information or to place an order call:

PHONE: +1 661-854-3166

FAX: +1 661-854-3850

email: customercare@salterlabs.com

salterlabs.com