

Below is information on our Rad-G. There are 2 models: a base model and one with integrated infrared temperature. Both have continuous and spot check modes, are super durable, rechargeable, measure through motion (wiggly kids), cold fingers and toes/low perfusion, and across all skin tones.


Both can come with either a reusable mini clip or a cable for adhesive sensors. There is also a cable that can be purchased separate if you'd like to have the option of both reusable mini clip and adhesives.

As for adhesive sensors, we have both the fold over and wrap around styles. **PLEASE SEE THE WEIGHT RANGES IN THE SENSOR DESCRIPTION**


RAD-G w/ TEMP, w/ ADHESIVE CABLE (PRODUCT #9895)

9895	<p>Rad-G® with Temperature (with Patient Cable) Includes SpO2, Pulse Rate, Perfusion Index (Pi), pleth variability index (PVi), Respiration Rate from the Pleth (RRp), Temperature and one RD SET G15 patient cable. Sensor must be ordered separately. <i>Available in countries where approved.</i></p>	EA	
------	---	----	--

RAD-G w/TEMP, REUSABLE MINICLIP (PRODUCT #9210)

PN	Description	UM	Product Image
9210	<p>Rad-G with Temperature</p> <p>(with Sensor) Includes SpO₂, Pulse Rate, Perfusion Index (Pi), pleth variability index (PVi), Respiration Rate from the Pleth (RRp), Temperature and one Rad-G reusable sensor.</p> <p><i>Available in countries where approved.</i></p>	EA	


RAD-G w/ REUSABLE MINI CLIP, BASE MODEL (PRODUCT #9847)

9847	<p>Rad-G Pulse Oximeter (with Sensor)</p> <p>Includes SpO₂, Pulse Rate, Perfusion Index (Pi), pleth variability index (PVi), Respiration Rate from the Pleth (RRp) and one Rad-G reusable sensor.</p>	EA	
------	---	----	--

RAD-G w/ ADHESIVE CABLE-BASE MODEL (PRODUCT #9849)

<p>9849</p>	<p>Rad-G Pulse Oximeter (with Patient Cable) Includes SpO₂, Pulse Rate, Perfusion Index (Pi), pleth variability index (PVi), Respiration Rate from the Pleth (RRp) and one RD SET G15 patient cable. Sensor must be ordered separately.</p>	<p>EA</p>	 <p>The image shows a white Rad-G Pulse Oximeter device with a black screen displaying vital signs: 97% SpO₂, 74 bpm Pulse Rate, and a pleth waveform. To the right is a coiled white patient cable with a connector at one end and a sensor at the other.</p>
-------------	---	-----------	--

ACCESSORY CABLE-Purchase in addition to either model that comes with the reusable miniclip (PRODUCT #4773)

<p>4773</p>	<p>RD SET G15-05 RD SET Patient Cable, 5 ft For use with Rad-G devices and RD SET series sensors</p>	<p>EA</p>	 <p>The image shows a coiled white patient cable with a connector at one end and a reusable miniclip sensor at the other.</p>
-------------	---	-----------	---

YI CABLE-DIRECT CONNECT OR SENSOR CONNECTION




4054

RD SET YI


Multi-site Reusable Sensor, 3 ft.
1/box, non-sterile
Multiple foam and adhesive wraps
Weight >1 kg
Not made with natural rubber latex.

EA




4653	<p>Rad-G YI Sensor Multisite Reusable Sensor, 3 ft. For direct connection with Rad-G Pulse Oximeter 1/box, non-sterile Weight > 1 kg Includes 2 foam wraps Not made with natural rubber latex</p>	EA	
1602	<p>Foam Wrap For use with YI Sensors 12 / pack, non-sterile Not made with natural rubber latex</p>	Pack	
1600	<p>Accessory Adhesive Squares For use with YI Sensors 144 / pack, non-sterile</p>	Pack	


NEO ADHESIVE SENSORS (ALSO USED FOR ADULTS)

4003	<p>RD SET Neo Adult / Neonatal SpO₂ Adhesive Sensor Single-Patient-Use 20 / box, non-sterile Weight < 3 kg or > 40 kg Not made with natural rubber latex Use Replacement tape 4674 Parameters Supported: SpO₂, PR, Pi, RRp, PVi</p>	Box	 <p>The image shows the RD SET Neo SpO₂ Adhesive Sensor. It consists of a white cable with a connector at one end and a tan adhesive sensor pad at the other. The sensor pad is rectangular with rounded corners and features a small black square sensor window. The word 'Neo' is printed vertically on the pad, along with a starburst symbol at the top and the Masimo logo at the bottom.</p>
------	---	-----	--


INFANT ADHESIVE SENSORS

4002	<p>RD SET Inf Infant SpO₂ Adhesive Sensor Single-Patient-Use 20 / box, non-sterile Weight 3 kg – 20 kg Not made with natural rubber latex Use Replacement tape 4046 Parameters Supported: SpO₂, PR, Pi, RRp, PVi</p>	Box	 <p>The image shows the RD SET Inf SpO₂ Adhesive Sensor. It consists of a white cable with a connector at one end and a tan adhesive sensor pad at the other. The sensor pad is rectangular with rounded corners and features a small black square sensor window. The word 'Inf' is printed vertically on the pad, along with a starburst symbol at the top and the Masimo logo at the bottom.</p>
------	--	-----	---

PEDIATRIC ADHESIVE SENSORS

4001	<p>RD SET Pdt</p> <p>Pediatric SpO₂ Adhesive Sensor Single-Patient-Use 20 / box, non-sterile Weight 10 kg – 50 kg Not made with natural rubber latex Parameters Supported: SpO₂, PR, Pi, RRp, PVi</p>	Box	
------	---	-----	---

ADULT ADHESIVE SENSORS

4000	<p>RD SET Adt</p> <p>Adult SpO₂ Adhesive Sensor Single-Patient-Use 20 / box, non-sterile Weight > 30 kg Not made with natural rubber latex Parameters Supported: SpO₂, PR, Pi, RRp, PVi</p>	Box	
------	--	-----	---